ENGLISH

Addendum to the ReNew Laparoscopic Hand Pieces Instructions for Use for Single Patient Use Disposable Tips

INDICATIONS

The ReNew single patient use disposable tips are to be used with the ReNew Laparoscopic Hand Pieces and they are indicated for cutting, grasping, dissecting and coagulation of tissue in endoscopic and laparoscopic surgical procedures.

CONTRAINDICATIONS

The ReNew single patient use disposable tips are not intended for use except as indicated.

DEVICE DESCRIPTION

Tips are available in various configurations and they are provided sterile, as single patient use, disposable.

INSTRUCTIONS FOR USE

Insert the threaded end of the disposable tip into the distal end of the shaft of the hand piece and hand tighten the tip clockwise until it is tight and secure.

WARNINGS AND CAUTIONS

- A thorough understanding of the principles and techniques involved in laparoscopic and electrosurgical procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments.
- The product is sterile unless the package is opened or damaged.
- Do not use if package is opened or damaged. Examine this device prior to use, do not use if damage is found.
- Do not reuse or resterilize.
- This single patient use device has not been designed with the intent of reprocessing (cleaning, disinfecting/sterilizing) and used on another patient
- Reprocessing and/or resterilization may compromise the integrity of the product which may result in its malfunctioning and unintended patient injury
- The physical characteristics and quality of this device may be adversely affected and the device may not remain safe and effective for its intended use when reprocessed and/or resterilized

Following are the risks associated with reuse:

- Risk of mechanical failure from undue stresses
- Risk that cleaning and decontamination may be ineffective, increasing the likelihood of cross contamination or infection
- Risk of residue from decontamination agents
- Risk of material alteration due to the decontamination process, including the use of chemical agents and elevated temperatures

WARRANTY INFORMATION

Microline Surgical warrants that its instruments are free from any defects in both material and workmanship. Microline Surgical shall not be held liable for any incidental or consequential damage of any kind. Work performed on an instrument will void this warranty. Gross abuse or neglect of a Microline Surgical instrument will void this warranty.

This product is an original MICROLINETM product, manufactured for single patient use only by Microline Surgical, Inc. Microline Surgical, Inc. is the owner of, and has exclusive rights to, the MICROLINE trademark. Unauthorized use of single patient use MICROLINE products, as well as any false designation of origin, or false or misleading description or representation of fact on or in connection with the MICROLINE trademark, will be considered to constitute trademark infringement and/or unfair competition.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician

Medical Device Safety Service GmbH
Schiffgraben 41
D-30175 Hannover
Germany

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