Microline Surgical, Inc.
50 Dunham Road, Suite 1500, Beverly, MA 01915 U.S.A.
Tel. +1(978) 922-9810 | Fax: +1(978)922-9209
Web: www.microlineurgical.com | Email: info@microlineurgical.com

09-39-00124-00-ART Rev J 01/2019

ENGLISH

ReNew Reusable Laparoscopic Instruments Tips

Instructions For Use

INDICATIONS FOR USE
ReNew Reusable Laparoscopic Instrument Tips are intended to be used with ReNew Reusable Laparoscopic Instrument Handpieces, and are indicated for cutting, grasping, dissecting, and coagulation of tissue in endoscopic and laparoscopic surgical procedures.

CONTRAINDICATIONS

None known.

DEVICE DESCRIPTION

The **ReNew Reusable Laparoscopic Instrument Tip** is intended to be individually attached to ReNew laparoscopic handpieces and used to manipulate tissue by mechanical means during use in surgical procedures. The device may be connected via a ReNew handpiece to a legally marketed electrocautery unit/radio frequency (RF) generator for dissection and coagulation of tissue. The device may be connected via a ReNew handpiece to a U.S. FDA-cleared electrocautery generator for dissection and coagulation of tissue. The device is detachable and may be interchangeably installed on all models of ReNew 5mm diameter handpieces. The device is intended to be reusable after processing/re-processing (cleaning, decontamination, packaging and sterilization). Consult the ReNew handpiece instructions for use for further information.

- Not made with natural rubber latex.
- Not made with bisphenol-A.
- Not made with Vinyl or Poly Vinyl Chloride (PVC).
- Not made with phthalates.

HOW SUPPLIED

ReNew Reusable Laparoscopic Instrument Tips are supplied non-sterile and must be sterilized before use.

WARNINGS

- Failure to adequately process and sterilize the device could lead to infection.
- The device is not intended to be bent, used to pry, or for any other use not approved by the U.S. FDA.
- Bending or using the device to pry could cause the device to fail and/or harm the patient.
- No modification of this device is allowed.
- Not intended for use with silicone or oil-based lubricants.
- Not intended for use with non-Microline, third party handpieces.
- A thorough understanding of the principles and techniques involved in endoscopic/laparoscopic procedures is essential to avoid harm to the patient or operator and to ensure the device is used as intended.
- The device has been validated for 25 sterilization cycles. Actual device performance may vary depending on care, use and handling of the device.
- Do not use device if package is opened or damaged. Examine the device prior to use. Do not use if damage is found.

WARNINGS WHEN ACTIVATING FOR ELECTROSURGERY

- Be aware that all exposed metal on the device is capable of treating tissue when the device is activated. Use caution to avoid inadvertent treatment of tissue.
- A thorough understanding of the principles and techniques involved in electrological procedures is essential to avoid shock and burn hazards to both patient and medical personnel or to avoid damage to the device or other medical instruments.
- Skin-to-skin contact (for example between the arms and body of the patient) should be avoided, for example by insertion of dry gauze.
- The patient should not come into contact with metal parts which are earthed or which have an appreciable capacitance to earth (for example operating table supports, etc.).
- Use the device with caution in the presence of pacemakers, as electrological devices may cause interference with pacemakers or other active implants.
- Place any monitoring electrodes being used as far away as possible from the device to avoid electrical interference with monitoring equipment.
- Avoid need-to-monitor electrodes.
- Use monitoring systems incorporating high frequency current limiting devices.
- DO NOT use electrocurey in the presence of flammable anesthetics or other flammable gases, near flammable fluids or objects, or in the presence of oxidizing agents, as fire could result.
- The device is designed for use in a dry area. The device is not intended to be used in a humid area, as increased humidity may result in injury or electrical shock.
- Insulation failures may result in burns or other injuries. Visual inspection alone may not be sufficient to confirm that the insulation is intact, and electric strength testing should be additionally considered.
- The exposed metal surface of the device may remain hot enough to cause burns after the radio frequency (RF) current is deactivated.
- Install device on the ReNew handpiece only when not connected to the electrological unit (ESU) or when the energy is off. Failure to do so may result in an injury or electrical shock.
- Aspirate fluid pools from the area being activated prior to activating RF current to the device. Conductive fluids (e.g., blood or saline) in direct contact with or in close proximity to the device may carry electrical current or heat away from tissue, which may cause unintended burns.
- Use the lowest setting possible on the RF Generator to achieve the desired tissue effect to avoid overtreatment, which could result in swelling, fluid, seroma or unintended tissue necrosis.
- Use only with a legally-marked, U.S. FDA cleared ESU.

PRECAUTIONS

- Follow all instructions, warnings, and precautions provided with the device and all associated devices.
- The device must be removed from the ReNew handpiece before cleaning and sterilization.
- Assembly and disassembly of the device with ReNew handpieces is intended to be achieved by hand without tools. Use of tools can damage the electrical insulation on the device and on the shaft of the handpiece.
- Any deviations from the processing instructions and parameters for cleaning, sterilization and reprocessing of the device as provided herein must be validated by the user.

ASSEMBLY AND DISASSEMBLY

- Manually attach the ReNew handpiece handles and the device jaws into their fully closed position.
- Insert the threaded end of the device into the distal end of the ReNew handpiece shaft.
- With one hand, hold the insulated rear-hub portion of the device in place and with the other hand turn the handpiece rotation knob clockwise until the rear hub of the tip is screwed tight and flush up against the shaft of the hand piece until the insulated portions of the device with ReNew handpieces is intended to be achieved by hand without tools.
- Rotate the ReNew handpiece handle to verify proper opening and closing of the device.
- To disassemble begin by squeezing the handle to fully close the jaws of the device.
- With one hand, hold the insulated rear-hub portion of the device in place and with the other hand turn the handpiece rotation knob counter clockwise until the rear hub of the device is fully unscrewed and separated from the shaft of the hand piece.
- For additional assistance and clarification refer to the illustrations and directions given in the instructions for use for the ReNew reusable handpiece.

MICROLINE HAS VALIDATED THE FOLLOWING CLEANING AND STEAM STERILIZATION PARAMETERS FOR USE WITH ReNew REUSABLE LAPAROSCOPIC INSTRUMENT TIPS

PREPROCESSING INSTRUCTIONS

- The device must be removed from the handpiece before cleaning and sterilization.
- MANUAL CLEANING (minimum required steps)**
 - Rinse with running tap water while actuating moving parts to remove gross soil.
 - Prepare enzymatic detergent solution in accordance to the manufacturer’s recommendations.
 - Immerse in enzymatic solution, actuate all moving parts, allow to soak for a minimum of five (5) minutes.
 - Clean the device of all visible gross contamination using soft brushes.
 - Rinse under running lukewarm tap water while actuating all moving parts.
 - Rinse under running deionized water while actuating all moving parts. Repeat 3 times.
 - Wipe with a clean, soft, lint-free cloth. Visually inspect to ensure all soil is removed.
- Note: Prior to sterilization, a high-temperature water-soluble lubricant must be applied to all movable parts of the device.
- Silicon- or oil-based lubricants are not recommended for use as these have not been validated and may prevent sterilization.

AUTOMATED CLEANING (minimum required steps)

- Rinse with running tap water while actuating moving parts to remove gross soil. If necessary, manually clean the device of any visible signs of contamination using a soft brush. Load device(s) into automated washer. Load device(s) in the fully open position, and, if a tray is used, do not cover the tray.
- Automated wash cycle parameters as follows:
Instrument Cycle: High

Phase	Time (Minutes)	Temperature	Detergent
Pre-Wash	2:00	Cold Tap Water	N/A
Wash 1	2:00	65 °C	Alkaline Detergent
Rinse 1	00:15	Hot Tap Water	N/A
PURW Rinse	1:00	90 °C	N/A
	0:00	98.8 °C	N/A
- Visually inspect to ensure all soil is removed.
- Package appropriately for sterilization.
- In addition to proper and thorough cleaning (see instructions) make sure the device is properly lubricated between cases:

 - If a lubricant is not added during the automated wash, it is important prior to sterilization that a water-soluble lubricant is added to aid and ensure a smooth articulation for all the movable parts of the device.

German Manual Pre-Cleaning:

- * Note that the Manual Pre-Cleaning method documented herein is intended to reflect a minimum requirement, this should be conducted before the automated cleaning method.
 - Rinse the tips under running tap water while actuating moving parts to remove gross soil.
 - Brush tips with a soft-bristled brush to aid in the removal of soil during the rinse.
 - Prepare detergent bath of Alkaline detergent per manufacturer’s instructions in a sonication unit.
 - Actuate the tips in the detergent bath before the sonication is started.
 - Sonicate the tips in the detergent bath for 10 minutes.
 - Place the tips into the tray or into the open position, and place the lid on the tray.
 - Transfer the tips onto the appropriate rack system contained inside the washer for processing.

German Automated Cleaning:

- Reference VARIO TD program as applicable.
- Note that the Automated cleaning method is intended to reflect a minimum requirement:
 - Select the appropriate cycle as listed below:
Motor speed: HIGH

PHASE	TIME (MINUTES)	TEMPERATURE	DETERGENT	
Pre-wash 1	02:00	Cold tap water	N/A	
Wash 1	05:00	Drain	55° C	Alkaline
Rinse 1	03:00	Drain	Tap water at <50° C	N/A
Final Rinse (PURW rinse)	02:00	Drain	RODI/Tap water <50° C	N/A
 - Disinfect: “The A₀ concept” (ref. DIN EN ISO 15883) and its national requirements must be taken into consideration in the on-site process validation, carried out by the user.
 - Visually inspect to ensure all soil is removed.

STERILIZATION PARAMETERS

Sterilize with the device in the open position. Device may be sterilized in a U.S. FDA cleared sterilization tray.

Sterilization Method	Temperature	Time	Dry Time	Sterilization Wrap
Pre-vacuum	132 °C	4 minutes	30 minutes	Wrapped in two layers of 1-ply polypropylene wrap of appropriate basis weight to accommodate the weight of the device or tray, using sequential envelope folding techniques.
Gravity	132 °C	15 minutes	20 minutes	

SPECIFIC STERILIZATION PARAMETERS FOR EUROPEAN COUNTRIES

Microline Surgical recommends the following parameters as minimum sterilization cycle parameters for Europe, except for France and Switzerland:

Sterilization Cycle	Exposure Time at 132 °C	Dry Time	Sterilization Wrap
Pre-vacuum	3-5 minutes	30 minutes	Wrapped in two layers of 1-ply polypropylene wrap of appropriate basis weight to accommodate the weight of the device or tray, using sequential envelope folding techniques.

For France and Switzerland, the following minimum sterilization cycle parameters are recommended:

Sterilization Cycle	Exposure Time at 132 °C	Dry Time	Sterilization Wrap
Pre-vacuum	18 minutes	30 minutes	Wrapped in two layers of 1-ply polypropylene wrap of appropriate basis weight to accommodate the weight of the device or tray, using sequential envelope folding techniques.

STORAGE CONDITIONS

Pre-sterilization: Store in a dry location.

Post-sterilization: Store per user specific protocols for storage of terminally sterilized devices.

REPROCESSING INSTRUCTIONS

- When a hand or instrument is removed from the ReNew handpiece. (See disassembly instructions above)
- Immerse device in sterile distilled water.
- When a device is being transferred to Central Sterile Processing, it should be wrapped in a towel moistened with distilled water.
- Thoroughly clean the device according to instructions provided in the Cleaning section above.
- Sterilize the device according to instructions provided in the Sterilization section above.
- Conduct a visual and functional inspection of the device per the Assembly and Disassembly instructions and Warnings above.
- If any component appears to be damaged, do not use the device.

DISPOSAL OF DEVICE

To reduce contamination risks, the device shall be disposed of in compliance with all applicable local laws and regulations. CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

CUSTOMER SERVICE CONTACT INFORMATION

For more information or assistance, please contact Microline Surgical at the address, phone, or email listed above.

WARRANTY

Microline Surgical warrants that its devices are free from any defects in both material and workmanship at the time of shipment. Microline Surgical shall not be held liable for any incidental or consequential damage of any kind or from damage resulting from neglect, abuse or mishandling of the device.
Work performed on a device will void this warranty.

Use with non-Microline, third-party handpieces will void this warranty.

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