Addendum to the ReNew Laparoscopic Hand Pieces Instructions for Use for Single Patient Use Disposable Tips**INDICATIONS**

The ReNew single patient use disposable tips are to be used with the ReNew Laparoscopic Hand Pieces and they are indicated for cutting, grasping, dissecting and coagulation of tissue in endoscopic and laparoscopic surgical procedures.

CONTRAINDICATIONS

The ReNew single patient use disposable tips are not intended for use except as indicated.

DEVICE DESCRIPTION

Tips are available in various configurations and they are provided sterile, as single patient use, disposable.

INSTRUCTIONS FOR USE

Insert the threaded end of the disposable tip into the distal end of the shaft of the hand piece and hand tighten the tip clockwise until it is tight and secure.

WARNINGS AND CAUTIONS

- A thorough understanding of the principles and techniques involved in laparoscopic and electrosurgical procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments.
- The product is sterile unless the package is opened or damaged.
- Do not use if package is opened or damaged. Examine this device prior to use, do not use if damage is found.
- Do not reuse or resterilize.
- This single patient use device has not been designed with the intent of reprocessing (cleaning, disinfecting/sterilizing) and used on another patient
- Reprocessing and/or resterilization may compromise the integrity of the product which may result in its malfunctioning and unintended patient injury
- The physical characteristics and quality of this device may be adversely affected and the device may not remain safe and effective for its intended use when reprocessed and/or resterilized

Following are the risks associated with reuse:

- Risk of mechanical failure from undue stresses
- Risk that cleaning and decontamination may be ineffective, increasing the likelihood of cross contamination or infection
- Risk of residue from decontamination agents
- Risk of material alteration due to the decontamination process, including the use of chemical agents and elevated temperatures

RETURN POLICY

Returns must be made within 30 days of shipment and must be in the original, unopened packaging. Disposable products must be returned in complete, unopened cases. Returned product is subject to a 25% restocking charge based on the original purchase price. Prior to returning any device, contact Microline Surgical Customer Service to obtain a Returned Goods Authorization Number (RGA). Contact Microline's Customer Service department or visit (www.microlinesurgical.com) for information. Unauthorized returns will not be accepted. All returns must be shipped freight prepaid. The following conditions are to be met when returning a used single use device:

- Proper cleaning and sterilization if used or a declaration of decontamination
- The device(s) must be assigned an RGA number
- The RGA number must be clearly visible on the outside of the box it is shipped in

WARRANTY INFORMATION

Microline Surgical's single patient use tips are unconditionally guaranteed against defects in material and workmanship. Microline shall not be held liable for any incidental or consequential damage of any kind. This warranty is valid only to the original purchaser of the device for the period of one (1) year. Microline will, at its option and without charge, replace tips which Microline determines to be defective in material or workmanship when used for their intended surgical purpose.

Microline Surgical, Inc. does not warrant reconditioned or reconstructed single patient use medical devices, in full or in part, that have been or will be resold, modified or treated by any third party, and accepts no responsibility or liability for such reconditioned or reconstructed devices.

NOTE:

This product is an original MICROLINE™ product, manufactured for single patient use only by Microline Surgical, Inc. Microline Surgical, Inc. is the owner of, and has exclusive rights to, the MICROLINE trademark. Unauthorized use of single patient use MICROLINE products, as well as any false designation of origin, or false or misleading description or representation of fact on or in connection with the MICROLINE trademark, will be considered to constitute trademark infringement and/or unfair competition.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician

CONTACT INFORMATION

Microline Surgical, Inc.
800 Cummings Center, Suite 166T
800 Cummings Center, Suite 157X
Beverly, MA 01915, USA
Tel: 978-922-9810
Fax: 978-922-9209

Medical Device Safety Service GmbH
Schiffgraben 41
D-30175 Hannover, Germany

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