M/L-10 Multi-Fire Clip Applier Cartridge Instructions For Use

Microline Surgical, Inc.

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50 Dunham Road, Suite 1500

Beverly, MA 01915 U.S.A.

Tel: +1(978) 922-9810

Fax: +1(978) 922-9209

Web: www.microlinesurgical.com

Email: info@microlinesurgical.com

09-39-00443-06-ART Rev B Nov, 2023

ENGLISH

M/L-10 Multi-Fire Clip Applier Cartridge – Instructions For Use

INDICATIONS FOR USE

The M/L-10 Clip Applier is indicated for occluding and ligating vessels, ducts, tracts and other tubular structures during laparoscopic and general surgical procedures.

Notes:

- Refer to M/L-10 Multi-Fire Clip Applier Instructions For Use, 09-39-00444-00-ART for additional instructions.
- Additional IFUs may be obtained free of charge by contacting Microline Customer service, or by visiting <https://www.microlinesurgical.com/resources/ifus/>

CONTRAINDICATIONS

- The M/L-10 Clip Applier is not intended for use except as indicated.
- The titanium clips are not intended to be used in the heart or blood vessels that are part of the central circulatory system.

WARNING

- If during the inspection of the shipping carton, any of the above items show any signs of damage, do not use the instrument.
- A thorough understanding of the principles and techniques involved in laparoscopic and electrosurgical procedures is essential to avoid shock and burn hazards to both patient and medical personnel and damage to the device or other medical instruments.
- The product is sterile unless the package is opened or damaged.
- Do not use if package is opened or damaged. Examine this device prior to use, do not use if damage is found.
- Do not reuse or resterilize.
- Not intended for use with non-Microline, third party devices.
- This single use device has not been designed with the intent of reprocessing
- Reprocessing and/or resterilization may compromise the integrity of the product which may result in its malfunctioning and unintended patient injury
- The physical characteristics and quality of this device may be adversely affected and the device may not remain safe and effective for its intended use when reprocessed and/or resterilized
- If an adverse event occurs using this device please report it to Microline Surgical Customer Service as well as the Competent Authority in your Member State (EU).

MR Compatibility

Non-clinical testing has demonstrated that the M/L-10 clips are MR Conditional. These devices can be scanned safely under the following conditions:

- Static magnetic field of 1.5-Tesla (1.5T) or 3.0-Tesla (3.0T).
- Spatial gradient field of up to:
 - 6,020 G/cm (60.20 T/m) for 1.5T systems
 - 3,010 G/cm (30.10 T/m) for 3.0T systems.
- Maximum whole body averaged specific absorption rate (SAR) of:
 - 4.0 W/kg for 15 minutes of scanning at 1.5T.
 - 4.0 W/kg for 15 minutes of scanning at 3.0T.

1.5T RF heating

In non-clinical testing with body coil excitation, the ligation clips produced a temperature rise of 1.0°C or less at a maximum whole body averaged specific absorption rate (SAR) of 4.0 W/kg, as assessed by calorimetry for 15 minutes of scanning in a 1.5T Siemens Espree (MRC30732) MR scanner with SYNGO MR B17 software.

3.0T RF heating

In non-clinical testing with body coil excitation, the ligation clips produced a temperature rise of less than 1.0°C or less at a maximum whole body averaged specific absorption rate (SAR) of 4.0 W/kg, as assessed by calorimetry for 15 minutes of scanning in a 3.0T Siemens Trio (MRC20587) MR scanner with SYNGO MR A30 4VA30A software.

Caution: The RF heating behavior does not scale with static field strength. Devices which do not exhibit detectable heating at one field strength may exhibit high values of localized heating at another field strength.

MR Artifact

In testing using a 3.0T system with gradient-echo sequencing, the shape of the image artefact produced follows the approximate contour of the device and extends radially up to 9 mm from the implant.

PRE-CHECK INSTRUCTIONS

Carefully inspect the packaging to ensure that there are no signs of damage. If damage is visible, do not use the device. The device is:

- Not made with natural rubber latex.
- Not made with bisphenol-A.
- Not made with Vinyl or Poly Vinyl Chloride (PVC).
- Not made with phthalates.

OPERATION NOTE

Before firing a clip, verify that the ligation site is free of any obstructions.

1. Verify that there is a clip properly loaded in the jaws prior to inserting the clip applier through the cannula.
2. Position the clip carefully around the vessel to be ligated.
3. With full visualization of the ligation site, squeeze the trigger all the way so that the trigger touches the handle to close the clip around the vessel. The trigger will not reopen until the trigger is beyond the 7th click (see step 5 below).
4. Releasing the trigger will open the clip applier jaws and automatically engage the next sequential new clip. The clip applier is now ready for re-firing.
5. For:
 - a. the cholangiogram procedure, whereby a clip is used to temporarily secure the catheter inside the common bile duct without occluding the catheter, the degree of clip closure can be adjusted. This is accomplished by releasing the trigger any time after the 7th click. This earlier

release capability allows for the desired amount of clip closure.

- b. maximum vessel occlusion (clip fully closed), simply depress instrument trigger fully so that the trigger and handle are touching and then fully release so that the subsequent clip can be engaged in the jaws.

CAUTIONS

1. If all the clips are to be used, be advised that the last two clips are colored blue to warn the surgeon that the cartridge is almost out of clips. The clip applicator should be removed from the cavity once the last colored (blue) clip has been fired. The cartridge should then be removed and discarded.
2. Engaging the jaws over a previously applied clip should be avoided, as it may cause the handpiece to misfire.
3. The application of clips over a staple line has not been validated for use with M/L-10 devices.
4. The cartridge must be removed from the handpiece after use, or prior to the cleaning and sterilizing processes. Microline Surgical shall not be responsible for any damage occurring because of non-compliance with this warning.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE (SSCP)

The Summary of Safety and Clinical Performance (SSCP) is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI. Eudamed can be accessed at <https://ec.europa.eu/tools/eudamed>. The Basic UDI-DI for the M/L-10 Clip Cartridges is 081109901DOC002S7

ENVIRONMENTAL PROTECTION

To reduce contamination risks, the device with unused clips shall be disposed of in compliance with all applicable local, state and federal laws and regulations.

Cleaning, sterilization, and grinding to enable recycling may be considered environmentally friendly, and supportive of circularity. Alternatively, Incineration may also be considered a safe disposal method.




WARRANTY INFORMATION

Microline Surgical warrants that its instruments are free from any defects in both material and workmanship. Microline Surgical shall not be held liable for any incidental or consequential damage of any kind. Work performed on an instrument will void this warranty. Gross abuse or neglect of a Microline Surgical instrument will void this warranty.

This product is an original MICROLINE™ product, manufactured for single patient use only by Microline Surgical, Inc. Microline Surgical, Inc. is the owner of, and has exclusive rights to, the MICROLINE trademark. Unauthorized use of single patient use MICROLINE products, as well as any false designation of origin, or false or misleading description or representation of fact on or in connection with the MICROLINE trademark, will be considered to constitute trademark infringement and/or unfair competition.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

SYMBOL DEFINITIONS

	Manufacturer
	Catalog Number
	Lot Number
	Use By Date
	Do not use if package is damaged
	Global Trade Item Number
	Sterilized by Ethylene Oxide
	Date of Manufacturer / Made in USA
	MR Conditional
	Authorized Representative for Switzerland
	Sterile Barrier / Method of Sterilization
	Consult instructions for use
	Authorized representative in the European Community
	Conformity European. This symbol means that the device conforms with European directives and regulations for product health, safety, and environmental protection.
	Do Not Resterilize
	Do not re-use
	Quantity
Rx Only	Prescription Only
	Not made with natural rubber latex.
	Medical Device
	UK Conformity Assessed
	Authorized representative in the United Kingdom
	For additional language IFUs scan this QR code.

Medical Device Safety Service GmbH
Schiffgraben 41
D-30175 Hannover
Germany

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MDSS CH GmbH