ENTceps / MicroCeps / CardioForceps Instructions For Use

Microline Surgical, Inc. 迈克兰器械有限公司.
50 Dunham Road, Suite 1500
Beverly, MA 01915 U.S.A.
Tel: +1(978) 922-9810 | Fax: +1(978) 922-9209
Web: www.microlinesurgical.com
Email: info@microlinesurgical.com

09-39-00450-00-ART, Rev. C 04/2021

ENGLISH

ENTceps / MicroCeps / CardioForceps – Instructions For Use

For use only with Microline PowerPack™ Dual Control Footswitch (REF 309-004R) and Microline PowerPack™ UPS (Universal Power Supply) (REF 200-003R) or Microline UPS (Universal Power Supply) (REF 200-004R, 200-006R)

DEVICE DESCRIPTION

The ENTceps / MicroCeps / CardioForceps are a footswitch-activated, hand-held surgical instrument that is powered by a Universal Power Supply (UPS). They contain a heating element at the distal tip of the instrument. The heating element is activated by means of the reusable PowerPack Dual Control Footswitch. The power cable extends from the handle of the instrument and terminates at the Instrument connector for electrical connection to the UPS.

INTENDED USE

The devices are single use and are intended to be used once only for a single patient. The devices are intended for the simultaneous cutting and cauterization of soft tissue during surgery. The devices may also be used for cutting natural or synthetic, non-metallic, sutures during surgery.

1. Remove the device from its sterile packaging.

Note: Follow the appropriate Instructions For Use for the PowerPack Dual Control Footswitch and UPS used for the procedure.

2. Insert the connector from the PowerPack Dual Control Footswitch (non-sterile) into the Footswitch connector receptacle of the UPS unit as per Footswitch Instructions For Use. Ensure the connector plug key properly into the receptacle. Uncoil the power cord, pass the connector end of the power cord off the sterile field. Insert the connector of the power cord firmly into the instrument connector receptacle of the UPS. Ensure connector plugs key properly into receptacle. The device is now ready for use.

3. Pre-check the Microline system for correct set-up and operation. Set-up the UPS and PowerPack Dual Control Footswitch following the applicable Instructions for Use. Place the device's tips in a small amount of sterile saline. Close the device and activate the device on high power by use of the PowerPack Dual Control Footswitch. There should immediately be an audible hissing sound as the saline is heated by the element, and the UPS will emit a high pitch audible signal tone. Repeat with the PowerPack Dual Control Footswitch operated on variable and the control knob on the UPS turned to the number 8 position. A lower pitch audible signal tone will be emitted. Repeat pre-check by attempting to activate the device by closing the forceps without stepping on the foot pedal. If there is a hissing sound as the saline is heated by the element, do not use the Microline system. Contact your Microline representative for repair or replacement.


Note: The heating element and tone will not operate unless the device is fully closed, the PowerPack Footswitch is connected to the UPS, and a pedal is depressed.

Note: The footswitch variable power output is adjustable when the device and PowerPack Footswitch are connected to the UPS and the footswitch is activated with the left footpedal. (See UPS Instructions for Use)

4. To seal/cut tissue, grasp the desired tissue between the tips of the device. Squeeze the handle grips fully together to apply pressure to the tissue. Depress either pedal of the footswitch. Depressing the left footpedal will activate the heating element at the tip of the device on variable power. (See Note above) An audible low-frequency tone will accompany the activation of the heating element. Depressing the right pedal on the footswitch will activate the heating element of the tip of the device on high power. An audible high-frequency tone will accompany the activation of the heating element. Generally, lower heat ranges increase the sealing capabilities and increase the time to cut tissue. Higher heating ranges decrease the time to cut and may decrease seal integrity.
5. After the desired sealing/cutting of the tissue, open the device and then release the pedal of the footswitch. This deactivates the heating element.
Note: After removing the device, examine tissue for hemostasis. If hemostasis is not present, use appropriate techniques to achieve hemostasis.
6. If desired, progress to a new region of tissue to be sealed/cut.
Note: It may be desirable to clean the tips of the device during the surgical procedure. A saline-moistened surgical gauze or sponge may be applied gently along length of the tips of the forceps to remove buildup of coagulated blood and tissue debris. Avoid twisting or pulling on tips. Applying excess pressure can bend the tips or cause tip misalignment.
7. At the end of the surgical procedure, disconnect and discard the device.


ENTceps (sterile)

- 1 – Instrument Connector
- 2 – Tip Cleaning Direction
- 3 – Heating Element


MicroCeps (sterile)

- 1 – Instrument Connector
- 2 – Tip Cleaning Direction
- 3 – Heating Element
- 4 – Finger Location (Squeeze Here)


CardioForceps (sterile)

- 1 – Instrument Connector
- 2 – Tip Cleaning Direction
- 3 – Heating Element

Dual Control Footswitch Non Sterile

- 1 – Variable Power Pedal (left)
- 2 – High Power Pedal (right)
- 3 – Power Supply Connector

Universal Power Supply Non Sterile

Label graphics may vary

PRECAUTIONS AND WARNINGS

- The product is sterile unless the package is opened or damaged.
- Do not use if package is opened or damaged. Examine this device prior to use, do not use if damage is found.
- The devices are not designed for reuse or reuse decontamination processes. Reuse of single-use devices creates a potential risk of patient or user infections, injury, illness or death.
- Do not use if the instrument or the cord is damaged.
- Do not use a scalpel or other sharp metal instrument to clean the device. Doing so may damage the tips and could prevent the device from functioning properly.
- Refrain from unnecessary activation of the heating element while there is no tissue grasped within the device. This activity may result in premature degradation of the device. Device is not intended for continuous use. A typical duty cycle of approximately five (5) seconds on, ten (10) seconds off is recommended.
- Do not immerse the device's handle in liquids.
- Do not touch electrosurgical (Bovie) electrode to the device.
- Use the device only with the Microline PowerPack Dual Control Footswitch and UPS Unit. Use of any other power supply may damage the device and could prevent it from functioning properly during use.

- Procedures using instruments for sealing/cutting should be performed only by persons having adequate training and familiarity with these surgical techniques. Consult the medical literature relative to techniques, complications and hazards prior to performance of any procedure. Surgeons using this device should be familiar with the specific anatomy of the region in which they intend to perform the procedure.
- There are no unusual risks associated with the proper disposal of this equipment. It shall be disposed of in compliance with all applicable local, state and federal laws and regulations
- There are no user-serviceable parts.
- If an adverse event occurs using this device please report it to Microline Surgical Customer Service as well as the Competent Authority in your Member State (EU).
- Does not contain latex.

Warning:

- Do not use in the presence of flammable materials (e.g. alcohol, flammable anesthetics).
- Always disconnect the instrument before discarding; the PowerPack Dual Control Footswitch and UPS are reusable.

CONTRAINDICATIONS

The ENTceps / MicroCeps / CardioForceps are not to be used as a fallopian tube sterilization device.

COMPLIANCE WITH STANDARDS

When used with the UPS, the device complies with IEC60601-1 requirements for type BF applied part and meets electromagnetic compatibility requirements of IEC60601-1-2.

HANDLING, STORAGE AND CLEANING

Do not drop. Store in a cool, dry place. Avoid prolonged exposure to extreme temperatures. Do not immerse the device in liquids.

LIMITED WARRANTY

Microline Surgical warrants that its instruments are free from any defects in both material and workmanship. Microline Surgical shall not be held liable for any incidental or consequential damage of any kind.

Work performed on an instrument will void this warranty.

Gross abuse or neglect of a Microline Surgical instrument will void this warranty.

SYMBOL DEFINITIONS

	Manufacturer		Consult instructions for use
	Catalog Number		Authorized representative in the European Community
	Batch code		Conformity European. This symbol means that the device conforms with European directives and regulations for product health, safety, and environmental protection.
	Temperature Limit		Prescription Only
	Do not use if package is damaged		Sterilized by Gamma Radiation
	Global Trade Item Number	Made in USA of US and imported parts.	
	Date of Manufacture		Atmospheric Pressure Limitation
	Variable Power		Use By Date
	Humidity Limitation		Do not re-use

Medical Device Safety Service GmbH
Schiffgraben 41
D-30175 Hannover
Germany

FDA cleared K000250

Manufactured under one or more of US patents: 7,033,351; 6,908,463; 6,860,880; 6,695,837; 6,626,901. Patents Pending. All Rights Reserved.

Recyclable and printed on paper containing post-consumer recycled fiber.

Australian Sponsor:

Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

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