M/L-10
MULTI-FIRE
CLIP
APPLIER

INSTRUCTIONS FOR USE

TABLE OF CONTENTS

ENGLISH4

FRANÇAIS9

ESPAÑOL14

ITALIANO19

DEUTSCH24

SIMPLIFIED CHINESE29

CHINESE33

M/L-10 REUSABLE MULTI- FIRE CLIP APPLIER

INSTRUCTIONS FOR USE

The M/L-10 Multi-Fire Clip Applier consists of:

1. One reusable autoclaveable handpiece – CAT. #1001/1002.
2. One removable "Stainless Steel" outer tube – CAT. #1201(included with CAT. #1001/1002).
3. Optional sterilization tray(s) - CAT. #3705 (**sold separately**) & 3707 (included w/clip applier & sold separately).
4. One sterile disposable back-loading slide cartridge - (**The cartridge, containing Titanium medium large clips, is sold separately.**)

WARNING:

*If during the inspection of the shipping carton, one of the above items shows any sign of damage **DO NOT USE.** Please contact your local distributor.*

LOADING CLIP APPLIER

1. Peel the clip cartridge pouch open and remove it from the pouch.
2. Fully squeeze the handpiece trigger against the handle. (Fig. 1)
3. Align the arrow on the back of the clip applier with the arrow of the cartridge as shown in fig 3 (Arrow of the cartridge should be seen on the "non-clip" side of the cartridge).
4. Insert the cartridge all the way in to the back of the handpiece. (Fig. 2), making sure it is handled by its sides or proximal end behind "the ladder" in order to avoid any displacement of clips or ladder.
5. Release the trigger. (This will immediately engage a clip into the handpiece jaws).
6. The instrument is ready to be used.


Fig. 1


Fig. 2


Fig. 3

OPERATION

NOTE:

Before firing a clip, verify that the ligation site is free of any obstructions.

1. Verify that there is a clip properly loaded in the jaws prior to inserting the clip applier through the cannula.
2. Position the clip carefully around the vessel to be ligated.

3. With full visualization of the ligation site, squeeze the trigger all the way so that the trigger touches the handle to close the clip around the vessel. The trigger will not reopen until the trigger is beyond the 7th click (see 5 below).
4. Releasing the trigger will open the clip applicator jaws and will automatically engage the next sequential new clip. Now the clip applicator is ready for re-firing.
5.
 - a. For the cholangiogram procedure, whereby a clip is used to temporarily secure the catheter inside the common cystic duct without occluding the catheter, the degree of clip closure can be adjusted. This is accomplished by releasing the trigger anytime after the 7th click at the ratchet mechanism. This earlier release capability allows for the desired amount of clip closure.
 - b. For maximum vessel occlusion (clip fully closed) simply depress instrument trigger fully so that the trigger and handle are touching and then fully release so that the subsequent clip can be engaged in the jaws.

CAUTIONS:

- A. If all the clips are to be used, the last **two colored clips will warn the surgeon that cartridge is almost out of clips**. The clip applicator should be removed from the cavity once the last colored clip has been fired. The cartridge should then be removed and disposed of.
- B. Engaging the jaws over a previously applied clip should be avoided as it may cause the handpiece to misfire

WARNING:

Make sure cartridge is removed from the handpiece after use, or prior to cleaning and sterilizing process. Microline PENTAX shall not be responsible of any damage occurring as a result of a non-compliance with this warning.

DISASSEMBLING THE CLIP APPLIER HANDPIECE

To disassemble clip applicator for reprocessing follow the steps below:

1. Remove and discard the disposable clip cartridge by fully squeezing the trigger against the handle as shown on Fig.4. If all 19 clips have not been fired squeezing the handle and removing the cartridge will engage and release 1 clip. In this case, make sure to fire the engaged clip onto something such as a drape so that the contaminated clip does not fall onto the floor, stick to the jaws or fall inside the instrument shaft potentially causing a jam or safety hazard. *(This step serves as a warning note on the "Cleaning and Sterilization" steps. Microline PENTAX shall not be responsible of any damages occurring as a result of a non-compliance with this warning.*
2. Remove the outer tube (Fig. 6) by fully squeezing the trigger and unscrewing the knob clockwise (Fig. 5).


Fig. 4


Fig. 5


Fig. 6

Notes: Use only Microline PENTAX's sterile disposable back-loading slide cartridge. Microline PENTAX will not be responsible for any damage caused by other clips.

Micro line PENTAX™

Advanced Surgical Instrumentation

RECOMMENDED CLEANING & STERILIZATION STEPS:

Warning:

Prior to cleaning the clip applicator instrument, make sure that the disposable cartridge is removed and discarded. Microline PENTAX shall not be responsible for any damages occurring as a result of a non-compliance with this warning.

Note:

The M/L-10 applicator should be stored and reprocessed in a protective sterilization tray.

1. Keeping the instrument moist:

Immediately after surgical use, remove the shaft's outer tube. Immerse instruments in a tray, basin or appropriate container of sterile distilled water. When the clip applicator handpiece and its outer tube are being transferred to Central Sterile Processing, they should be wrapped in a moistened towel of distilled water.

2. Enzymatic soaking & ultrasonic wash:

The instrument and its outer tube should be immersed in an ultrasonic bath filled with a mixture containing 60mL of ENZOL® Enzymatic Detergent per 4 L of water @ 22-40°C. This will aid in the cleaning, flushing and removal of any debris. Sonicate device for 10 minutes.

3. Rinsing

After the soaking cycle, the clip applicator and its outer tube should be removed from the enzyme solution and rinsed with regular tap water.

4. Cleaning

Once thoroughly rinsed, the clip applicator and its outer tube should be immersed in an approved cleaning solution. The handle, shaft and rotation knob should be cleaned with a soft cloth or sponge soaked in concentrated ENZOL® Enzymatic Detergent until no further bioload is observed on brush or towel. The distal end of the instrument must be cleaned with a "soft bristle" brush to avoid any damage.

5. Rinsing

Once the washing cycle is completed, the clip applicator handpiece and its outer tube should be rinsed in normal tap water, and then rinsed 3 times with deionized water. Wipe instrument with a clean soft cloth.

6. Magnifying Lens Inspection

After the hospital cleaning procedure is completed, an inspection of the instrument should be done to insure it is thoroughly cleaned.

7. Drying

Prior to wrapping instruments for storage, they must be completely dry. Any residual moisture could result in oxidation and corrosion that could affect the performance, function, and useful life of the instrument.

8. Lubricating

It is important, prior to sterilization, that a water-soluble lubricant is employed to preserve the articulation and action of all movable parts of the instrument.

Re-Assembling

9. Re-assemble the outer tube on the instrument. To perform this, squeeze the handpiece trigger, mount the outer tube over the clip applicator shaft, and then secure it by screwing the tube knob counterclockwise until fully threaded.

Note:

Your reprocessing protocol may require that the applicator and the outer tube be reassembled after the sterilization process.

Micro line PENTAX™

Advanced Surgical Instrumentation

STERILIZATION:

Warning:

Prior to sterilizing the clip applier instrument, make sure cartridge is removed. Microline PENTAX shall not be responsible for any damages occurring as a result of a non-compliance with this warning.

1. Prior to sterilization, the instrument must be thoroughly cleaned.
2. Wrap the instrument or use sterilization tray (Cat #3705/3707).
3. Sterilize the instrument following the protocol provided by the sterilizer manufacturer.

Microline PENTAX recommends the following as minimum sterilization cycle requirements:

Gravity Cycle:	3 minutes @ 270°F (132°C)
Pre-Vacuum Cycle:	3 minutes @ 270°F (132°C)

INSPECTION:

NOTE:

M/L-10 Inventory Inspection Procedure (to be performed once a month)

An abbreviated inspection should be performed in the OR each time before handing the M/L-10 to the surgeon for use:

- 1. Check jaws for misalignment, damage and caked on bio- burden.
- 2. Check shaft for straightness. Remove outer tube and check inner shaft for damage and caked on bio-burden.
- 3. Squeeze and release handle 20 plus times to check for sticking. Check handle for damage. Never pry handle and trigger open.
- 4. Fire a minimum of 2 clips to check clip closure, engagement into jaws and clip/jaw alignment.
- 5. In addition to proper and thorough cleaning (see instructions) make sure the M/L-10 is properly lubricated between cases:
 - » A. Fully immerse the instrument in an approved instrument lubricant.
 - » B. If unable to fully immerse the instrument make sure to thoroughly lubricate both jaw and rotation knob area.
 - » C. After lubrication process perform step 3 and 4 of the inspection if the applier handle was sticking or locking.

WARRANTY:

Microline PENTAX warrants that this instrument is free from defects in both material and workmanship.

Microline PENTAX shall not be held liable for incidental or consequential damage of any kind. This warranty is valid only to the original purchaser of the instrument and for a period of one (1) year.

Work performed on this instrument by anyone other than a Microline PENTAX Authorized Service Center will void this warranty.

SERVICE AND REPAIR:

Prior to returning this instrument for repair, call your Microline PENTAX Distributor to obtain authorization.

The instrument will be returned un-repaired to the sender if the following conditions are not met:

- The instrument must be given a "Return Goods Authorization" (R.G.A) number.
- The instrument must be cleaned and sterilized prior to returning for repair.
- The R.G.A. Number must be clearly visible on the outside of the box it is shipped in.

**Micro
line PENTAX™****Advanced Surgical Instrumentation**

For any questions, please call your local Microline PENTAX distributor or Microline PENTAX Customer Service at:

CONTACT INFORMATION:**Microline PENTAX, Inc.**

800 Cummings Center

Suite 157X

Beverly, MA 01915

U.S.A.

Tel: (978) 922-9810

Fax: (978) 922-9209

Web: www.microlinepentax.comE-mail: info@microlinepentax.com**EU REPRESENTATIVE:****PENTAX Europe GmbH**

Julius-Vosseler-Str. 104

22527 Hamburg • Germany

Tel: +49-40-56 192 0

Fax: +49-40-56 042 13

IMPORTANT MESSAGE TO OUR VALUED CUSTOMERS

Congratulations on your purchase of the Microline PENTAX Reusable M/L -10 Multi-Fire Clip Applier. This precision, reusable instrument performs a vital surgical function and should be handled with the utmost care. **We strongly encourage you to use the protective Microline PENTAX M/L-10 Sterilization Tray, Catalog No. 3707.**

For questions related to the care, maintenance, and handling of the M/L-10, please contact your authorized Microline PENTAX Representative. If you are unsure as to your Representative, contact our website **www.microlinepentax.com** or Microline PENTAX's main office at: +1 (978) 922-9810.

THANK YOU FOR CHOOSING MICROLINE PENTAX

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