MICROLINE™ SURGICAL

ReNew Electrocautery Probe Handpiece and Electrode Tips Instructions For Use

Microline Surgical, Inc. 迈克兰器械有限公司
50 Dunham Road, Suite 1500
Beverly, MA 01915, U.S.A.
电话: (978) 922-9810
传真: (978) 922-9209
网站: www.microlinesurgical.com
电子邮件: info@microlinesurgical.com

09-39-45-ART Rev K 01/2019

ENGLISH

ReNew Electrocautery Probe Handpiece and Electrode Tips

INDICATIONS

The ReNew Electrocautery Probe and electrode tips are intended to cauterize tissues in endoscopic and laparoscopic surgical procedures.

CONTRAINDICATIONS

The ReNew Electrocautery Probe and electrode tips are not intended for use except as indicated.

INSPECTION

Carefully inspect the shipping carton and its contents for any sign of damage. If damage is visible, DO NOT USE the instrument.

PRIOR TO USE

This instrument is delivered without electrode tips. Electrode tips are available in various configurations.

NOTE: The ReNew electrode tips are specifically designed to be used only with the ReNew Electrocautery Probe.

The handpiece is delivered with a detachable cap that allows for connection of the generator cord and also flushing of the lumen, when removed from the handpiece. Ensure that the cap is securely tightened to the handle prior to use.

HANDPIECE/TIP ASSEMBLY

1. Insert an electrode tip into the distal end of the shaft, and screw the tip in clockwise until tip is screwed tight. A tip is fully assembled to the shaft when the white marker band located on the distal end of the shaft is no longer visible.
2. The instrument is ready for use.

PRECAUTIONS

1. CAUTION: Prior to use, assure that the plastic hub of the electrode tip is fully in contact with the insulation tube of the handpiece and that no gap exists between the two parts at that contact point.
2. Do not use the handpiece if the threaded distal end of the shaft appears worn, damaged or missing.
3. Do not use mechanical means to assemble/disassemble the electrode tip to/from the device. Serious malfunction of the device may occur.
4. Use only with safety certified generators.

ELECTROCAUTERY

First, connect the electrosurgical cord (not furnished with this instrument) to the handpiece by placing the female end of the cord on the male adapter pin on the handpiece, which is recessed in the cap. Plug the other end of the cord into the monopolar receptacle of the generator.

PRECAUTIONS

1. A complete understanding of the principle of monopolar electrocautery surgical procedures is necessary to avoid accidental shocks, burns, or potential gas embolism to the patient.
2. Keep the power as low as possible to achieve the desirable effect.
3. Do not activate the generator until the tip is in contact with the tissue or is in a position to deliver high frequency energy to the tissue.
4. Be sure that a return electrode grounding pad has been properly attached to the patient and the generator.
5. The ReNew Electrocautery Probe has an integrated sealing feature at the distal end of the shaft, which is an integral part of the system.
6. Electrical Ratings: To reduce risk of electrical burn or shock, hand pieces and tips are not to be used beyond their maximum recurring peak voltage ratings: 1KVp in cutting mode, 2.35KVp in coagulation burst mode, 4.0 KVp in spray coagulation mode.

HANDPIECE CLEANING

CAUTION: ELECTRODE TIPS MUST BE REMOVED FROM THE ELECTROSURGICAL CAUTERY PROBE HANDPIECE PRIOR TO CLEANING AND STERILIZATION OF THE HANDPIECE SINCE THE TIPS ARE SINGLE USE DISPOSABLE.

Manual Cleaning:

Note: For cleaning purposes, a device consists of shaft and cautery cap.

1. Disassemble device:
 - a. Remove the cautery cap from the proximal end of the device.
 - b. Ultrasonic cleaning:
 - i. Place device into ultrasonic bath filled with enzymatic detergent solution prepared in accordance with manufacturer's recommendations.
 - ii. Sonicate device for 10 minutes.
 - iii. Rinse device with running tap water to remove gross soil for one (1) minute.
2. Presoak in enzymatic Solution:
 - a. Prepare enzymatic detergent solution in accordance to the manufacturer's recommendations.
 - i. Flush any lumen-like channels with detergent, by inserting a 60mL luer type syringe filled with prepared detergent solution where the cap was removed and flushing with a minimum of 60 mL of detergent solution. Repeat filling the syringe and flushing the shaft for a total of three (3) times and 180 mL.
 - ii. Immerse device in enzymatic solution; allow to soak for a minimum of one (1) minute.

- b. After the soak time, while still immersed, use a soft bristled brush to remove all visible soil from the device. Pay attention to any challenge areas while brushing for a minimum of one (1) minute.
 - c. After the soak time, while still immersed, use an appropriately sized lumen brush to brush the interior of all lumens for a minimum of one (1) minute.
3. Lukewarm tap water rinse:
- a. Remove the device and rinse under lukewarm running tap water for a minimum of one (1) minute.
 - i. Insert a 60mL luer type syringe filled with lukewarm tap water where the cap was removed and flush with a minimum of 60 mL. Repeat filling the syringe and flushing the shaft for a total of three (3) times and 180mL. Ensure that clear water exits the other side while flushing. (CAUTION: Water temperatures >45°C may cause scalding).
4. Reverse Osmosis/ deionized RO/DI Rinse:
- a. Rinse the device in Reverse Osmosis/ deionized (RO/DI) water. Repeat two (2) additional times for a total of three (3) RO/DI rinses.
 - b. Insert a 60mL luer type (or equivalent) syringe filled with warm deionized (RO/DI) water where the cap was removed and rinse with 60 mL of deionized water. Repeat filling the syringe and flushing the shaft three (3) times for a total of 180mL.
 - c. Dry the device completely using clean lint-free cloth.
 - d. Visually inspect the device in a well-lit area to ensure all surfaces are clean. Repeat cleaning steps if visible soil is still observed.
 - e. Inspect the device for functionality and package appropriately for sterilization.

Automated Cleaning:

1. Unthread any tip that is on the handpiece.
2. Remove the cautery cap from the proximal end of the device.
3. Rinse with running tap water for 30 seconds to remove gross soil.
4. Flush the device with lukewarm tap water while fully immersed in a bath of lukewarm tap water (CAUTION: Water temperatures >45°C may cause scalding).
5. Use a syringe to flush the devices with lukewarm tap water.
6. Use an appropriately sized lumen brush to clean all hard to reach areas.
7. Load device(s) into automated washer. Load device(s) disassembled.
8. Automated wash cycle parameters as follows:
9. Instrument Cycle: High

US:

Phase	Time (minutes)	Temperature	Detergent
Pre-Wash	2:00	Cold Tap Water	N/A
Wash 1	2:00	65.5°C	Alkaline Detergent
Rinse 1	00:15	Hot Tap Water	N/A
PURW Rinse	1:00	90°C	N/A
Dry	6:00	98.8°C	N/A

German Automated Cleaning:

- a. Reference VARIO TD program as applicable.
- b. Note that the Automated cleaning method is intended to reflect a minimum requirement:
- c. Select the appropriate cycle as listed below:
Motor speed: HIGH

PHASE	TIME (MINUTES)	TEMPERATURE	DETERGENT
Pre-wash 1	02:00	Cold tap water	N/A
Drain			
Wash 1	05:00	55°C	Alkaline
Drain			
Rinse 1	03:00	Tap water at <50°C	N/A
Drain			
Final Rinse (PURW rinse)	02:00	RO/DI water <50°C	N/A
Drain			

- d. Disinfect: "The A₀ concept (ref. DIN EN ISO 15883) and its national requirements must be taken into consideration in the on-site process validation, carried out by the user."
- e. Visually Inspect to ensure all soil is removed.

HANDPIECE STERILIZATION

1. Prior to sterilization, the handpiece must be thoroughly cleaned.
2. Wrap the handpiece with the cap removed from the shaft.

Microline SURGICAL recommends using the following validated steam sterilization cycles as guidelines:

Sterilization Cycle	Temperature	Exposure Time	Dry Time	Sterilization Wrap
Pre-vacuum	132°C	4 minutes	30 minutes	Wrapped in two layers of 1-ply polypropylene wrap of appropriate basis weight to accommodate the weight of the device, using sequential envelope folding techniques.
Pre-vacuum	134°C	3 -5 minutes	30 minutes	
Pre-vacuum	134°C	18 minutes	30 minutes	

WARRANTY

Microline Surgical warrants that its instruments are free from any defects in both material and workmanship. Microline Surgical shall not be held liable for any incidental or consequential damage of any kind. Work performed on an instrument will void this warranty. Gross abuse or neglect of a Microline Surgical instrument will void this warranty.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

EC REP **Medical Device Safety Service GmbH**
Schiffgraben 41
D-30175 Hannover
Germany

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