# MICROLINE<sup>TM</sup> SURGICAL

## ReNew Electrocautery Probe Handpiece and Electrode Tips Instructions For Use

Microline Surgical, Inc.  
50 Dunham Road  
Beverly, MA 01915, U.S.A.  
Tel: (978) 922-9810  
FAX: (978) 922-9209

Web: [www.microlinesurgical.com](http://www.microlinesurgical.com)  
Email: [info@microlinesurgical.com](mailto:info@microlinesurgical.com)

09-39-45 Rev G 02/2015

# MICROLINE<sup>TM</sup> SURGICAL

ENGLISH

### ReNew Electrocautery Probe Handpiece and Electrode Tips

**INDICATIONS**  
The ReNew Electrocautery Probe and electrode tips are intended to cauterize tissues in endoscopic and laparoscopic surgical procedures

### CONTRAINDICATIONS

The ReNew Electrocautery Probe and electrode tips are not intended for use except as indicated.

### INSPECTION

Carefully inspect the shipping carton and its contents for any sign of damage. If damage is visible, DO NOT USE the instrument.

### PRIOR TO USE

This instrument is delivered without electrode tips. Electrode tips are available in various configurations. NOTE: The ReNew electrode tips are specifically designed to be used only with the ReNew Electrocautery Probe. The handpiece is delivered with a detachable cap that allows for connection of the generator cord and also flushing of the lumen, when removed from the handpiece. Ensure that the cap is securely tightened to the handle prior to use.

### HANDPIECE/TIP ASSEMBLY

1. Insert an electrode tip into the distal end of the shaft, and screw the tip in clockwise until tip is screwed tight. A tip is fully assembled to the shaft when the white marker band located on the distal end of the shaft is no longer visible.
2. The instrument is ready for use.

### PRECAUTIONS

1. CAUTION: Prior to use, assure that the plastic hub of the electrode tip is fully in contact with the insulation tube of the handpiece and that no gap exists between the two parts at that contact point.
  2. Do not use the handpiece if the threaded distal end of the shaft appears worn, damaged or missing.
  3. Do not use mechanical means to assemble/disassemble the electrode tip to/from the device. Serious malfunction of the device may occur.
  4. Use only with safety certified generators.
- ELECTROCAUTERY**  
First, connect the electrocautery cord (not furnished with this instrument) to the handpiece by placing the female end of the cord on the male adapter pin on the handpiece, which is recessed in the cap. Plug the other end of the cord into the monopolar receptacle of the generator.

### PRECAUTIONS

1. A complete understanding of the principle of monopolar electrocautery surgical procedures is necessary to avoid accidental shocks, burns, or potential gas embolism to the patient.
2. Keep the power as low as possible to achieve the desirable effect.
3. Do not activate the generator until the tip is in contact with the tissue or is in a position to deliver high frequency energy to the tissue.
4. Be sure that a return electrode grounding pad has been properly attached to the patient and the generator.
5. The ReNew Electrocautery Probe has an integrated sealing feature at the distal end of the shaft, which is an integral part of the system.
6. Electrical Ratings: To reduce risk of electrical burn or shock, hand pieces and tips are not to be used beyond their maximum recurring peak voltage ratings: 1KvP in cutting mode, 2.35KvP in coagulation burst mode, 4.0 KvP in spray coagulation mode.

### HANDPIECE CLEANING

**CAUTION: ELECTRODE TIPS MUST BE REMOVED FROM THE ELECTROSURGICAL CAUTERY PROBE HANDPIECE PRIOR TO CLEANING AND STERILIZATION OF THE HANDPIECE SINCE THE TIPS ARE SINGLE USE DISPOSABLE.**

1. Remove the cap from the proximal end of the device.
2. Presoak one device (shaft and cap) at a time in a solution containing 60mL of ENZOL® Enzymatic Detergent, or equivalent, per 4 L of water at 22-40°C.
3. Rinse the device (shaft and cap) with warm tap water to remove any visible residue and debris for one (1) minute.
4. Clean the device (shaft and cap) of any visible signs of contamination using a soft brush or soft towel soaked in concentrated ENZOL® Enzymatic Detergent until no further bio-load is observed on the brush or towel. Ensure that movable assemblies are advanced to their full range of motion for careful cleaning of all surfaces.
5. Rinse the device (shaft and cap) with warm tap water.
6. Insert a 60 mL "luer" type syringe filled with warm tap water where the cap was removed.
7. Flush the inside of the shaft and tip clean.
8. Repeat steps 6 & 7 three (3) times.
9. Flush inside of the shaft three (3) times with deionized water for a final rinse.
10. Rinse the device (shaft and cap) three (3) times with deionized water for a final rinse.

Note:

Other equivalent neutral pH cleaning solutions that are suitable for the product may be used. Always follow the manufacturer's instructions regarding concentration of the cleaning solution, temperature and exposure time.

- Mechanical Cleaning/Decontaminating/Disinfecting**
- After performing manual cleaning, the instrument may be placed in a suitable washer/decontaminator/disinfector basket and processed through a standard instrument washer/decontaminator/disinfector reprocessing cycle.
  - Follow the instructions provided by the manufacturer of the machine.
  - Use a suitable neutral pH cleaning agent according to the manufacturer's instructions.
  - After completion of the mechanical cleaning/decontaminating/disinfecting cycle, inspect all surfaces, cavities, lumens and openings for any visible debris.
  - Carry out additional manual cleaning if necessary.

### HANDPIECE STERILIZATION

1. Prior to sterilization, the handpiece must be thoroughly cleaned.
  2. Wrap the handpiece with the cap removed from the shaft.
- Microline SURGICAL recommends using the following validated steam sterilization cycles as guidelines:
- Gravity Cycle: 4 minutes @ 270°F (132°C)  
Pre-Vacuum Cycle: 30 minutes @ 250°F (121°C)  
Pre-Vacuum Cycle: 4 minutes @ 270°F (132°C)

### STERILIZATION PARAMETERS FOR EUROPEAN COUNTRIES

Microline SURGICAL recommends the following parameters as minimum sterilization cycle parameters for Europe, except for France and Switzerland:  
Pre-Vacuum Cycle: 3-5 minutes at 134°C - 137°C  
For France and Switzerland, the following minimum sterilization cycle parameters are recommended: Pre-Vacuum Cycle: 18 minutes at 134°C - 137°C, 2x102 kPa

### WARRANTY

Microline Surgical warrants that its instruments are free from any defects in both material and workmanship. Microline Surgical shall not be held liable for any incidental or consequential damage of any kind. Work performed on an instrument will void this warranty. Gross abuse or neglect of a Microline Surgical instrument will void this warranty.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

### CONTACT INFORMATION

Microline Surgical, Inc.  
50 Dunham Road  
Beverly, MA 01915 U.S.A.  
Tel: (978) 922-9810  
Fax: (978) 922-9209  
Web: [www.microlinesurgical.com](http://www.microlinesurgical.com)  
Email: [info@microlinesurgical.com](mailto:info@microlinesurgical.com)

Medical Device Safety Service GmbH  
Schiffgraben 41  
D-30175 Hannover  
Germany

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