# MICROLINE<sup>TM</sup>

## SURGICAL

### Gemini Multi-Fire Medium/Large and Large Clip Appliers Instructions for Use

Microline Surgical, Inc.  
800 Cummings Center, Suite 166T  
800 Cummings Center, Suite 157X  
Beverly, MA 01915, USA  
TEL: (978) 922-9810  
FAX: (978) 922-9209  
Web: [www.microlinesurgical.com](http://www.microlinesurgical.com)  
Email: [info@microlinesurgical.com](mailto:info@microlinesurgical.com)

09-39-00056-00 Rev. C 04/2010

ENGLISH

### Gemini Multi-Fire Medium/Large and Large Clip Appliers

#### INDICATIONS

The Gemini Medium/Large and Large Clip Appliers are indicated for occluding and ligating vessels, ducts, tracts and other tubular structures during laparoscopic and general surgical procedures.

#### CONTRAINDICATIONS

The Gemini Medium/Large and Large Clip Appliers are not intended for use except as indicated.

#### INSTRUCTIONS FOR USE

The Gemini Medium/Large and Large Clip Appliers consist of:

1. One reusable autoclavable handpiece
2. One removable stainless steel outer tube
3. One sterile disposable back-loading slide cartridge, containing Titanium large clips (sold separately).

#### WARNING


If during the inspection of the shipping carton, one of the above items shows any sign of damage, do not use the instrument.

#### PRODUCT OFFERINGS

The Gemini Clip Appliers are offered in a Medium/Large and Large size. Verify cartridge size to the clip applicator size prior to opening the cartridges. Verify color coordination between the cartridge and the clip applicator. Blue for Medium/Large and Green for Large. Use color matched components only.

#### LOADING CLIP APPLIER

1. Verify color coordination between the cartridge and the clip applicator. Blue for Medium/Large and Green for Large. Use color matched components only.
2. Peel the clip cartridge pouch open and remove it from the pouch.
3. Fully squeeze the handpiece trigger against the handle (Fig. 1).
4. Align the arrow on the back of the clip applicator with the arrow of the cartridge as shown in Fig. 3 (arrow of the cartridge should be seen on the "non-clip" side of the cartridge).
5. Insert the cartridge all the way into the back of the handpiece (Fig. 2), making sure it is handled by its sides or proximal end behind "the ladder", in order to avoid any displacement of clips or ladder.
6. Twist the cartridge cap to lock in place.
7. Release the trigger. (This will immediately engage a clip into the handpiece jaws).
8. The instrument is ready to be used.

#### OPERATION NOTE

Before firing a clip, verify that the ligation site is free of any obstructions.

1. Verify that there is a clip properly loaded in the jaws prior to inserting the clip applier through the cannula.
2. Position the clip carefully around the vessel to be ligated.
3. With full visualization of the ligation site, squeeze the trigger completely so that the trigger has a firm feel and nearly touches the handle, to close the clip around the vessel. The trigger will not reopen until the clip has reached the maximum size for a cholangiogram.
4. Releasing the trigger will open the clip applier jaws and will automatically engage the next sequential new clip. Now the clip applier is ready for re-firing.
5.
  - a. For the cholangiogram procedure, whereby a clip is used to temporarily secure the catheter inside the common bile duct without occluding the catheter, the degree of clip closure can be adjusted. This is accomplished by releasing the trigger soon after a distinctively louder click is heard (8<sup>th</sup> click). This earlier release capability allows for the desired amount of clip closure.
  - b. For maximum vessel occlusion (clip fully closed) simply depress the instrument's trigger fully so that the trigger and handle are touching and then fully release so that the subsequent clip can be engaged in the jaws.

#### CAUTIONS

- a) If all the clips are to be used, the last two colored clips will warn the surgeon that the cartridge is almost out of clips. The clip applier should be removed from the cavity once the last colored clip has been fired. The cartridge should then be removed and disposed.
- b) If the device is cycled after the last clip is fired, the trigger will lock. Remove the cartridge to clear the device. Do not over-power the lock as this may damage the handpiece.
- c) Engaging the jaws over a previously applied clip should be avoided as it may cause damage to the handpiece or cause it to misfire.

#### WARNING

Make sure that the cartridge is removed from the handpiece after use, or prior to cleaning and sterilization process. Microline shall not be responsible of any damage occurring as a result of a non-compliance with this warning.

#### DISASSEMBLING THE CLIP APPLIER HANDPIECE

To disassemble the clip applier for reprocessing follow the steps below:

1. Remove and discard the disposable clip cartridge by fully squeezing the trigger against the handle as shown on Fig.4. If all clips have not been fired, squeezing the handle and removing the cartridge will engage and release one clip. In this case, make sure to fire the engaged clip onto something such as a drape, so that the contaminated clip does not fall onto the floor, stick to the jaws or fall inside the instrument shaft potentially causing a jam or safety hazard. (This step serves as a warning note on the "Cleaning and Sterilization" steps. Microline shall not be responsible for any damages occurring as a result of a non-compliance with this warning).
2. Remove the outer tube (Fig. 6) by fully squeezing the trigger and unscrewing the knob clockwise (Fig. 5).


Fig. 4


Fig. 5


Fig. 6

**Notes:** Use only a Microline Gemini sterile disposable cartridge. Microline will not be responsible for any damage caused by other cartridges or clips.

### RECOMMENDED CLEANING & STERILIZATION STEPS

#### Warning

Prior to cleaning the clip applicator instrument, make sure that the disposable cartridge is removed and discarded. Microline shall not be responsible for any damages occurring as a result of a non-compliance with this warning.

#### CLEANING

- 1) Immediately after surgical use, remove the shaft's outer tube. The instrument and its outer tube should be soaked in a neutral pH enzymatic solution containing 60mL of ENZOL® Enzymatic Detergent (or equivalent) per 4 L of water @ 22-40°C.
- 2) After the soaking cycle, the clip applicator and its outer tube should be removed from the enzyme solution and rinsed with regular tap water.
- 3) The handle, shaft and rotation knob should be cleaned with a soft cloth or sponge soaked in concentrated ENZOL® Enzymatic Detergent until no further bioburden is observed on cloth or sponge. The distal end of the instrument must be cleaned with a "soft bristle" brush to avoid any damage.
- 4) The clip applicator handpiece and its outer tube should be rinsed in normal tap water, and then rinsed 3 times with deionized water for a final rinse. Wipe instrument with a clean soft cloth.
- 5) After the hospital cleaning procedure is completed, an inspection of the instrument should be done to ensure it is thoroughly cleaned.
- 6) It is important, prior to sterilization, that a water-soluble lubricant is used to preserve the articulation and action of all movable parts of the instrument.
- 7) Re-assemble the outer tube on the instrument. To perform this, squeeze the handpiece trigger, mount the outer tube over the clip applicator shaft, and then secure it by screwing the tube knob counterclockwise until it is fully threaded.

The hospital's reprocessing protocol may require that the clip applicator and the outer tube be reassembled after the sterilization process. Prior to wrapping the instrument for storage, it must be completely dry. Any residual moisture could result in oxidation and corrosion that could affect its performance, function and useful life.

**Note:** Other equivalent neutral pH cleaning solutions that are suitable for the product may be used. Always follow the manufacturer's instructions regarding concentration of the cleaning solution, temperature and exposure time.

#### Mechanical Cleaning/Decontaminating/Disinfecting

- After performing manual cleaning, the instrument may be placed in a suitable washer/decontaminator/disinfector basket and processed through a standard instrument washer/decontaminator/disinfector reprocessing cycle.
- Follow the instructions provided by the manufacturer of the machine.
- Use a suitable neutral pH cleaning agent according to the manufacturer's instructions.
- After completion of the mechanical cleaning/decontaminating/disinfecting cycle, inspect all surfaces, cavities, lumens and openings for any visible debris.
- Carry out additional manual cleaning if necessary.

#### STERILIZATION

##### Warning

Prior to sterilizing the clip applicator instrument, make sure that the cartridge is removed. Microline shall not be responsible for any damages occurring as a result of a non-compliance with this warning. Prior to sterilization, the instrument must be thoroughly

cleaned. Microline recommends the following validated sterilization cycles as guidelines:

Sterilizer Type:	Pre -Vacuum
Minimum Temperature:	270°F (132°C)
Cycle Time:	4 minutes
Minimum Dry Time:	N/A
Configuration:	Unwrapped or in sterilization case

Sterilizer Type:	Gravity
Minimum Temperature:	270°F (132°C)
Cycle Time:	15 minutes
Minimum Dry Time:	N/A
Configuration:	Unwrapped or in sterilization case

#### **STERILIZATION PARAMETERS FOR EUROPEAN COUNTRIES**

Microline recommends the following parameters as minimum sterilization cycle parameters for Europe, except for France and Switzerland:

Sterilizer Type:	Pre -Vacuum
Temperature:	134°C - 137°C
Cycle Time:	4-5 minutes
Minimum Dry Time:	N/A
Configuration:	Unwrapped or in sterilization case

For France and Switzerland, the following minimum sterilization cycle parameters are recommended:

Sterilizer Type:	Pre -Vacuum
Temperature:	134°C - 137°C, 2x10 <sup>2</sup> KPa
Cycle Time:	18 minutes
Minimum Dry Time:	N/A
Configuration:	Unwrapped or in sterilization case

#### **INSPECTION**

Gemini clip applicator Inventory Inspection Procedure (to be performed once a month). An abbreviated inspection should be performed in the OR each time before handing the Gemini clip applicator to the surgeon for use:

1. Check jaws for misalignment, damage and caked on bioburden.
2. Check shaft for straightness. Remove outer tube and check inner shaft for damage and caked on bioburden.
3. Squeeze and release handle 20 plus times to check for sticking. Check handle for damage. Never pry handle and trigger open.
4. Fire a minimum of 2 clips to check clip closure, engagement into jaws and clip/jaw alignment.
5. In addition to proper and thorough cleaning (see instructions) make sure the Gemini clip applicator is properly lubricated between cases:
  - a. Fully immerse the instrument in an approved instrument lubricant.
  - b. If unable to fully immerse the instrument make sure to thoroughly lubricate both jaw and rotation knob area.
  - c. After lubrication process perform step 3 and 4 of the inspection if the applicator handle was sticking or locking.

#### **WARRANTY**

Microline warrants that this instrument is free from defects in both material and workmanship. Microline shall not be held liable for incidental or consequential damage of any kind. This warranty is valid only to the original purchaser of the instrument and for a period of one (1) year. Work performed on this instrument by anyone other than a Microline Authorized Service Center will void this warranty.

#### **SERVICE AND REPAIR**

Prior to returning this instrument for repair, call your Microline distributor to obtain authorization. The instrument will be returned unrepaired to the sender if the following conditions are not met:

- The instrument must be given a "Return Goods Authorization" (RGA) number.
- The instrument must be cleaned and sterilized prior to returning for repair.
- The RGA number must be clearly visible on the outside of the box it is shipped in.

**Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

#### **CONTACT INFORMATION** **Microline Surgical, Inc.**

800 Cummings Center, Suite 166T  
800 Cummings Center, Suite 157X  
Beverly, MA 01915, USA  
Tel: (978) 922-9810  
Fax: (978) 922-9209

Medical Device Safety Service GmbH  
Schiffgraben 41  
D-30175 Hannover  
Germany

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