MICROLINE™

SURGICAL

Microline Surgical, Inc. 迈克兰器械有限公司.

50 Dunham Road, Suite 1500

Beverly, MA 01915 U.S.A.

Tel: +1(978) 922-9810 | Fax: +1(978) 922-9209 Web: www.microlinesurgical.com

Email: info@microlinesurgical.com

ENGLISH

PowerPack™ Dual Control Footswitch – Instructions For Use

Note: Additional IFUs may be obtained free of charge by contacting Microline Customer service, or by visiting <https://www.microlinesurgical.com/resources/ifus/>.

For use only with

Microline Cautery Instruments and Microline Model # 200-006R Universal Power Supply (UPS)

Device Description

The PowerPack Dual Control Footswitch Model # 309-004R (Footswitch) is a non-sterile, reusable unit intended to be used only with Microline Cautery Instruments and the Microline UPS. The Footswitch has one outlet cable which is connected to the UPS. The Footswitch has two power settings, variable and high. The amount of power delivered in the variable mode is determined by the button setting on the power supply. Variable power is activated by depressing the left footpedal. High power is activated by pressing the right footpedal (see figure).

PowerPack Dual Control Footswitch Model # 309-004R (Footswitch)

Intended Use

The Footswitch is intended for use only with Microline Cautery Instruments and the UPS for the simultaneous cutting and cauterization of soft tissue during surgery.

Operating Sequence

1. Refer to the appropriate Instructions For Use for both the UPS and the Microline Cautery Instrument being used.
2. Connect grounded, hospital grade power cord to the UPS AC connector receptacle.

3. Plug the power cord into a grounded hospital grade receptacle.
4. Align the key portion of the Microline Cautery Instrument connector with the UPS connector receptacle and insert.
5. Plug the Footswitch connector into the Footswitch connector receptacle of the UPS. Ensure the connector plug keys properly into the receptacle on the UPS.
6. Place the Footswitch on a level surface near enough to the surgeon that it can be activated by the surgeon's foot.
7. Turn on the UPS.
8. UPS: When the instrument heating element is activated, the UPS emits a tone that will last until the heating element is deactivated. The amount of heat delivered by the instrument heater may be adjusted by pressing the "+" or "-" push pads on the right side of the UPS. The numerical setting is shown on the LED bar graph above the push pads. Increasing the LED's numerical value setting increases the amount of heat delivered by the instrument when using a single-heat range instrument or the variable range when using the Footswitch accessory or multi-heat range instrument. Generally, lower heat ranges increase the sealing capabilities and increase the time required to divide tissue. Higher heat ranges decrease the time to divide and may decrease seal integrity.
9. When using Footswitch activated devices, depressing the Footswitch pedal will activate the heating element only after closing cautery device. See device specific instructions for use.
10. Activate the cautery Instrument by first closing the device and then depressing the left or right foot pedal (variable or high power respectively) depending on the desired-clinical result. When the heating element of the Cautery Instrument is activated the UPS emits an audible tone that lasts until the instrument heating element is deactivated. A low pitch tone is emitted when variable power is delivered. A high pitch tone is emitted when high power is delivered
11. At the end of the surgical procedure:
 - Turn the UPS off
 - Disconnect the Microline Cautery Instrument and
 - Disconnect the Footswitch.

Recommended Environmental Conditions

| | Transport | Storage | Operating |
|---|-------------------|-------------------|-------------------|
|  | -15°C to 40° | -15°C to 35° | 15°C to 35° |
|  | 10 to 95 | 10 to 93 | 30 to 93 |
|  | 700hPa to 1060hPa | 700hPa to 1060hPa | 700hPa to 1060hPa |

Precautions and Warnings

- Use only with Microline Cautery Instruments and UPS. Use of any other instrument may damage the Footswitch and could prevent it from functioning properly during use.
- Do not drop the Footswitch.
- Do not immerse the Footswitch in liquids.
- Do not sterilize the Footswitch.
- If using multiple Microline Cautery Instruments, do not allow instruments to contact each other. Do not touch electrosurgical (monopolar) electrodes to Microline Cautery Instruments.
- The Footswitch is not intended for continuous use. A typical duty cycle of approximately five (5) seconds on, ten (10) seconds off is recommended. If the Footswitch is used continuously, or the instrument connector becomes shorted, the output of the UPS may be interrupted. This will be indicated by an absence of audible tone from the UPS. Attempt to determine the cause of the problem before proceeding. If problem

cannot be resolved or power is not restored after several minutes of cool-down, contact Microline Surgical for instructions regarding repair or replacement of the device.

- There are no user-serviceable parts.
- If an adverse event occurs using this device, please report it to Microline Surgical Customer Service as well as the Competent Authority in your Member State (EU).

Warning:

- Do not use in the presence of flammable materials (e.g. alcohol, flammable anesthetics).
- The Footswitch is reusable.
- Avoid use in the proximity of other equipment that may interfere with the proper operation of the UPS.
- The Footswitch is intended to be used only in conjunction with the 200-006R UPS.
- No modification of the Footswitch is allowed.

Compliance with Standards

When used with the UPS, the Footswitch complies with IEC60601-1 (CAN/CSA-C22.2 No. 601.1-M90) requirements for type Class I, Type BF equipment, Type BF applied parts, and meets electromagnetic compatibility requirements of IEC60601-1-2.

Handling, Storage and Cleaning

- Do not drop.
- Store in a cool, dry place.
- Avoid prolonged exposure to extreme temperatures.
- Exterior of the Footswitch may be cleaned using a soft cloth moistened with a solution of water and a mild detergent or disinfectant (such as Cidex® or Sporidicin®).
- Do not immerse the Footswitch in liquids.

ENVIRONMENTAL PROTECTION

To reduce contamination risks, the Footswitch shall be disposed of in compliance with all applicable local, state and federal laws and regulations.

Cleaning, sterilization, and grinding to enable recycling may be considered environmentally friendly, and supportive of circularity. Alternatively, Incineration may also be considered a safe disposal method.

Limited Warranty

Microline Surgical warrants that its instruments are free from any defects in both material and workmanship. Microline Surgical shall not be held liable for any incidental or consequential damage of any kind.


Work performed on an instrument will void this warranty.

Gross abuse or neglect of a Microline Surgical instrument will void this warranty.

CAUTION

Federal law (USA) restricts this device to sale by or on the order of a physician.

Symbol Definition

| | | | |
|---|---|---|--|
|  | Manufacturer |  | Consult instructions for use |
|  | Catalog number |  | Quantity |
|  | Batch code | Rx Only | Prescription Only |
|  | Do not use if package is damaged |  | Non sterile |
|  | Global Trade Item Number | Made in USA of US and imported parts. | |
|  | Medical Device |  | Date of Manufacture / Made in USA |
|  | Date of Manufacture |  | High Power |
|  | Variable Power | | On (Power Connected to the main.) |
|  | Serial Number |  | Atmospheric Pressure |
|  | Temperature |  | Type BF applied Part |
|  | Relative Humidity |  | Off (Power disconnected from the main) |
|  | General Warning |  | Footswitch Connector |
|  | Authorized representative in the European Community |  | Conformity European. This symbol means that the device conforms with European directives and regulations for product health, safety, and environment protection. |
|  | Not made with natural rubber latex. |  | Authorized Representative for the United Kingdom |
|  | Authorized Representative for Switzerland |  | UK Conformity Assessed |

Medical Device Safety Service GmbH
Schiffgraben 41
D-30175 Hannover
Germany

MDSS-UK RP Ltd.
6 Wilmslow Road, Rusholme,
Manchester M14 5TP
United Kingdom

MDSS CH GmbH

Laurenzvorstadt 61, 5000 Arau, Switzerland

Manufactured under one or more of US patents: 7,033,351; 6,908,463; 6,860,880; 6,695,837; 6,626,901. Patents Pending. All Rights Reserved.

Australian Sponsor:
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

Recyclable and printed on paper containing post-consumer recycled fiber.
Cidex is a registered trademark owned by ASP Global Manufacturing, GmbH.
Sporicidin is a registered trademark owned by Contec, Inc.

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Cidex, ASP Global Manufacturing, GmbH'in tescilli ticari markasıdır.
Sporicidin Contec, Inc.'in tescilli ticari markasıdır.