

# REVEEL

ENDOSCOPIC RETRACTOR



**Radically Rethinking  
Endoscopic Retraction**






# INTRODUCING REVEEL

- GENTLE
- STRONG
- FAST

**GENTLE.** Surface with a Coefficient of Frictions (COF) 171% greater than traditional metallic retractors<sup>1</sup>. Using REVEEL results in: **less retraction force**; reduced possibility of inadvertent crushing of the liver; **less slippage**; and **less frustration** for the surgeon.

**ACTIVE TRACTION.** Patented retraction surface that really “sticks” to the liver and stops it from sliding around during the procedure. The combination of ACTIVE RETRACTION and the unique HARP SHAPE provides **strong** and truly **gentle** retraction.

**PERCUTANEOUS.** There is no need for additional trocars when using REVEEL. With the same diameter as an industry standard 5mm cannula, REVEEL’s profiled tip facilitates easy, low-force insertion, while the internal seals help maintain pneumoperitoneum and prevent messy fluid leaks.



**EXCELLENT VISUALIZATION.** Shaped for an unobstructed view, REVEEL offers exceptional exposure of the operating site and is optimized for large, fatty livers. The **harp shape** and rigid **central i-beam** enables the liver to be well-mobilized and prevent the left lobe from falling into the field of view.

## ORDERING INFORMATION FOR REVEEL:



**RV093001** SINGLE USE (6 Per Box)  
Percutaneous - inserted via a 5mm incision  
(heart width 90mm, working length 310mm)



**RT083001** SINGLE USE (6 Per Box)  
Trocars Compatible\* - 5mm trocar compatible,  
inserted via a trocar or percutaneously. (heart  
width 80mm, working length 300mm)

\*Compatible with industry leading 5mm trocars

**ERGONOMICS.** Comfortable, intuitive, and easy to use handle. Pull-and-go actuation allows the retraction surface to be locked in the desired shape. One-touch **QUICK RELEASE** mechanism allows for safe, rapid deactivation.



Scan the QR code to view REVEEL



[www.microlinesurgical.com](http://www.microlinesurgical.com)

**US**

50 Dunham Road Suite 1500 | Beverly, MA 01915 USA  
Phone 978-922-9810 | Fax 978-922-9209  
[customerservice@microlinesurgical.com](mailto:customerservice@microlinesurgical.com)

**EMEA**

Storkstraat 17B | 3833 LB Leusden, The Netherlands  
Phone +31 (0) 33 76 000 11 | +31 (0) 33 76 000 14  
[customerservicembv@microlinesurgical.com](mailto:customerservicembv@microlinesurgical.com)

© 2023 Microline Surgical, Inc. Microline is a trademark and/or registered trademark of Microline Surgical, Inc.

**References:** <sup>1</sup> Chiu P, Leung B, Sy S, Chan C, Wong S, Poon C (2013) The Safety and Efficacy of a New Liver Retractor for the Performance of Laparoscopic Upper Gastrointestinal Surgery – A Preclinical Evaluation Study. Dept. of Surgery, The Chinese University of Hong Kong, SAGES Scientific Session, April 2013

P/N: 09-39-00544-00 Rev B

**Important information:** Prior to use, refer to Instructions for Use supplied with each device for indications, contraindications, adverse events, suggested procedures, warnings and precautions for full product information.

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a Physician.